



# MultiChrome™


Veterinary Microbiology Culture System



## Your Key To Success In Controlling Mastitis

Mastitis is not only a physical problem for the herd, but an economic problem for the dairy producer whose prices are affected by the presence of this infection. MultiChrome™, with its unique chromogenic bi-plate media not only confirms the presence of Mastitis causing bacteria, but will separate them into gram positive and gram negative families. It then further identifies individual types of bacteria most associated with Mastitis by forcing them to ignite into distinct colors.

**In house Color Specific Detection And Identification of Six Common Mastitis Bacteria... In Just 24 Hours after incubation!**



S. agalactiae

S. uberis

S. aureus

Klebsiella

Citrobacter

E. coli

**Now Available! MultiChrome™ Mastitis “Quad” Plates for detection of infectious vs. environmental Mastitis causing bacteria.**



*“Kacey Products Providing Better Animal Care”*