




# MultiChrome™



## Veterinary Micro Incubator

**Digital incubator with twelve assignable, built in timers.  
Precise 24 hour incubation that alerts you when completed!**


The Kacey® Micro Incubator is the perfect companion to the MultiChrome™ microbiology system with its adjustable controlled heat at 98f / 37c and Patented, unique use of individual timers. No more guessing. Just *Load... Press Timer...and Walk Away!*

### MultiChrome Micro Incubator:

- Small foot print (8"w x 13"h x12"d)
- Precise temperature control
- Digital read out
- Trays cradle all forms of dishes
- Clear door for observation
- Twelve individual timers
- Audible/Visible alerts
- Energy efficient, UL listed
- Made for MultiChrome™ culture plates, but suitable for Sensitivity testing as well.



Each timer indicates **"Green"** during the 24 hr timing event and then turns to **"Red"** with audible sound when ready to read.



Contact us or your local Kacey Diagnostics dealer for more information about the MultiChrome™ Micro Incubator.

Authorized Dealer:



**KACEY™**  
DIAGNOSTICS

1854 A Hendersonville Rd  
Asheville, NC 28803  
Ph: 1(828)685-3569  
Fax: 1(828)685-3571  
[www.kaceyinc.org](http://www.kaceyinc.org)

*"Kacey Products Providing Better Animal Care"*