




Poison Test Station

COMPLETE "IN HOUSE" TESTING OF THREE COMMON POISONS

**Anti-Freeze
Ethylene Glycol**

**Rodenticide
Anti-Coagulant** **Insecticide/Pesticide
Organo-Phosphate**



Rodenticide Anti-Coagulant Test

Specifications:

- One step Go/No Go result
- No instrument required
- Stat results in under 5 minutes
- Screens for Rodenticides
- Screens Super Rodenticides
- 12 month expiration

Anti-Freeze Ethylene Glycol

Specifications:

- Feline sensitivity @ 20mg/dl
- Canine sensitivity @ 50mg/dl
- Only 20ul of plasma needed
- Detects 30 min after ingestion
- Results in 8 minutes
- Two color comparison chart
- 14 months expiration

Insecticide / Pesticide Organo-Phosphate

Specifications:

- Feline sensitivity 75 U/dl
- Canine sensitivity 125 U/dl
- Only 20ul of plasma needed
- Results in 6 minutes
- Easy color comparison chart
- 14 months expiration
- Low unit cost per strip

The Kacey Poison Control Station for Quick, Accurate results in **ONE** location!

Authorized Dealer:



KACEYTM
DIAGNOSTICS

1854 A Hendersonville Rd
Asheville, NC 28803
Ph: 1(828)685-3569
Fax: 1(828)685-3571

www.KaceyDiagnostics.com

Kacey Inc. © 2010 All Rights Reserved.

"Kacey Products Providing Better Animal Care"