



Veti-Fecal™

Veterinary Fecal Flotation System

Complete Kit With Stabilized, Non-Crystallizing Solution!

The Veti-Fecal is a complete program for ova identification for both floatation and centrifuge use. The solution does not crystallize in suspension enhancing microscopic confirmation and identification of ova. The zinc sulfate concentrate is stable for 24 months and working solution will maintain effectiveness for up to 18 months covering 50 tests!

Veti-Fecal™ features:

- Stable specific gravity of 1.18
- Enough concentrate for 50 tests
- Concentrate stable for 24 months
- Mixed solution stable for 18 months
- No suspended crystals in solution
- Homogenous solution after reconstitution
- Suitable for both floatation and centrifuge use
- Ergonomic bottle for easier pouring



Veti-Fecal Kit includes testing devices, collection instruction sheets, laboratory preparation instructions and enough concentrate to perform fifty tests.

Contact us or your local Kacey Diagnostics dealer for more information about the Veti-Fecal system.

Authorized Dealer:



KACEY™
DIAGNOSTICS

1854 A Hendersonville Rd
Asheville, NC 28803
Ph: 1(828)685-3569
Fax: 1(828)685-3571
www.kaceyinc.org

“Kacey Products Providing Better Animal Care”