


MultiChrome™ Veterinary Microbiology Culture System




Color Specific Detection And Identification of Common Pathogenic Bacteria...In 24 Hours!

MultiChrome is formulated to isolate Gram Positive and Gram Negative organisms. Unlike older media, this chromogenic media is designed to ignite individual, distinct colors within the bacteria allowing isolation and differentiation of pathogens for a fraction of the cost of an outside lab!.

MultiChrome can test for:

- Outer Ear Infections
- Urinary Tract Infections
- Conjunctival Infections
- Genital Infections
- Wound/abscess Pathogens
- Upper Respiratory Infections
- Skin Infections
- Mastitis
- Septicemia
- MRSA specific MultiChrome test available soon!


Simply inoculate bi-plate, place in incubator for 24 hrs & identify



Only a few of the bacteria detected!

Additional sensitivity testing of confirmed infections can be either run in-house or sent to a lab. MultiChrome has a 6 month refrigeration shelf life. In addition, we offer the Kacey Diagnostic MultiChrome Incubator with a small footprint and big features such as individual 24 hour timers for each test that is run.

Contact us or your local Kacey Diagnostics dealer for more information about the MultiChrome culture system.



Authorized Dealer:



KACEY™
DIAGNOSTICS

1854 A Hendersonville Rd
Asheville, NC 28803
Ph: 1(828)685-3569
Fax: 1(828)685-3571
www.kaceyinc.org

“Kacey Products Providing Better Animal Care”