




# MultiChrome™

## Veterinary Micro Incubator



**Digital incubator with twelve assignable, built in timers.  
Precise 24 hour incubation that alerts you when completed!**


The Kacey® Micro Incubator is the perfect companion to the MultiChrome™ microbiology system with its adjustable controlled heat at 98f / 37c and Patented, unique use of individual timers. No more guessing. Just *Load... Press Timer...and Walk Away!*

### MultiChrome Micro Incubator:

- Small foot print (8"w x 13"h x12"d)
- Precise temperature control
- Digital read out
- Trays cradle all forms of dishes
- Clear door for observation
- Twelve individual timers
- Audible/Visible alerts
- Energy efficient, UL listed
- Made for MultiChrome™, but suitable for Dermatophyte and Sensitivity testing.



Each timer indicates "Green" during the 24 hr timing event and then turns to "Red" with audible sound when ready to read.



Contact us or your local Kacey Diagnostics dealer for more information about the MultiChrome™ Micro Incubator.

Authorized Dealer:



**KACEY™**  
DIAGNOSTICS

1854 A Hendersonville Rd  
Asheville, NC 28803  
Ph: 1(828)685-3569  
Fax: 1(828)685-3571  
[www.kaceyinc.org](http://www.kaceyinc.org)

*"Kacey Products Providing Better Animal Care"*